

## Pour Plate Method: Best Practices

### Best Practices for Pour Plate Method

- According to a consensus of views, it is recommended to keep the molten agar in the water bath for no more than 3 to 4 hours; do not pour it until it has cooled to <math><50^{\circ}\text{C}</math> ( preferably <math>44^{\circ}\text{C}</math> to <math>46^{\circ}\text{C}</math>).
- The agar should not be melted more than one time.
- If it is necessary to dilute the suspension, use phosphate buffer pH 7.2.
- To decrease the risk of contamination, pour plates in a laminar-airflow cabinet. To decrease the risk of contamination when pouring multiple plates, flame the mouth of the flask before moving on to the next plate.
- The FDA BAM recommends filling plates with 12 ml to 15 ml of agar. The USP recommends a fill of 15 ml to 20 ml of agar.
- Some microorganism species may recover better on spread plates than pour plates. If this is the case, it is recommended to verify counts with a spread plate.
- Incubate most bacterial species for 48 to 72 hours. Note: Incubate *Candida albicans* and *Aspergillus brasiliensis* for 3 to 5 days.
- Some microorganism colonies, such as *Pseudomonas aeruginosa*, are very small. Count the colonies with the aid of an illuminated colony counter or magnifying glass.
- If determining the number of CFU per ml use the formula below.  
$$\text{Number of CFU per ml} = \# \text{ of Colonies} \times \text{Dilution factor} \times \frac{1}{\text{Aliquot}}$$
- Recovery will be lower on selective agar. If selective agar is used, test non-selective agar in parallel using the same microorganism suspension. A higher CFU concentration for the selective agar may be necessary.
- The value obtained when performing tests may differ from the mean assay value. Note: Microbiologics uses non-selective Tryptic Soy Agar when testing most microorganisms.


### References

- International Standard for Organization. 2009-02-01. Second Edition: Microbiology of food and animal feeding stuffs – Guidelines on preparations and production of culture media. Part 1
- Maturin, L. and Peeler, J. T. 2001. FDA Bacteriological Analytical Manual. Chapter 3, Aerobic Plate Count.
- Standard Methods for the Examination of Dairy Products. 2004. 17th Edition
- Standard Methods for the Examination of Water and Wastewater. 2005. 9020, Quality Assurance/Quality Control.
- United States Pharmacopeia. 2013. <61>, Microbiological Examination of Nonsterile Products: Microbial Enumeration Tests.

## ILLUSTRATED INSTRUCTIONS

Instructions for performing the Pour Plate Method

1




Place your molten agar in a circulating water bath set at <math>< 50^{\circ}\text{C}</math>, preferably

2




Prepare microorganism suspension according to protocol.

3




Pipette inoculum onto the bottom of the petri dish. Test each suspension in duplicate.

4




In accordance with your standards, pour 12-20 ml of molten agar over the inoculum.

5




Thoroughly mix the organisms throughout the agar by gently swirling the petri dish by alternate rotation.

6




After the agar has hardened, invert the plates.

7



Incubate.

8



Count the colonies. Colonies under the agar surface will be small.